

深圳市疾病预防控制中心 卫生检测报告

编号：B2020018

受检（委托）单位：深圳市疾病预防控制中心 Shenzhen Center for Diseases Control and Prevention	地址：南山区龙苑路8号 Nanshan District, No. 8 Longyuan Road
样品类别：鼻咽拭子 Nasopharyngeal swab	采（送）检单号： /
样品状态：液态 Liquid	采样方式： /
样品来源：送样 Sample	检测日期：2020年12月7日
检测项目：新型冠状病毒抗原检测试剂盒（胶体金法） SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	检测环境条件：温度：24℃ 湿度：42%
检测方法：深圳市绿诗源生物技术有限公司（新型冠状病毒抗原胶体金免疫层析法试剂盒说明书） Shenzhen Lvshiyuan Biotechnology Co., Ltd (SARS-CoV-2 Antigen Rapid Test Kit Colloidal Gold Chromatographic Immunoassay)	
检测结果： Detection results 检测过程选取核酸荧光 PCR 方法确诊阳性样本 100 份以及阴性样本 210 份，阳性符合率为 98%，阴性符合率为 100%，总体准确率为 99.35%。与其他相关病毒的交叉实验证明，该检测试剂与其他相关病毒无交叉反应。 In this clinical trial, 100 positive samples and 210 negative samples (confirmed by the nucleic acid fluorescence PCR method) were tested. The positive coincidence rate was 98%, the negative coincidence rate was 100%, and the total accuracy rate was 99.35%. The results of cross-experiments with other related viruses have proved that the test kit has no cross-reaction with other related viruses. 见下面附件 See the attachment below	

检测者：

复核者：

签发者：



签发日期：2020年12月7日（盖章）

Information of samples (Nasopharyngeal swabs storing in 3ml Hank's solution)							Brand: GreenSpring® Name: SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)
Number	Sample type	Sample collection date	Date of symptom onset	Ct value	Results (+/-)		
1	Nasopharyngeal swabs	2020/3/26	2020/3/24	22	++++		
2	Nasopharyngeal swabs	2020/1/15	2020/1/8	24.5	++		
3	Nasopharyngeal swabs	2020/3/28	2020/3/24	23	+++		
4	Nasopharyngeal swabs	2020/2/3	2020/1/27	23	+		
5	Nasopharyngeal swabs	2020/2/6	2020/1/20	22	+++		
6	Nasopharyngeal swabs	2020/4/9	2020/4/1	22	++++		
7	Nasopharyngeal swabs	2020/1/25	2020/1/17	23	++		
8	Nasopharyngeal swabs	2020/1/29	2020/1/27	23	+++		
9	Nasopharyngeal swabs	2020/2/15	2020/2/9	22	++++		
10	Nasopharyngeal swabs	2020/2/14	2020/2/14	22	++++		
11	Nasopharyngeal swabs	2020/3/20	2020/3/17	28	+		

12	Nasopharyngeal swabs	2020/3/26	2020/3/26	28	±
13	Nasopharyngeal swabs	2020/3/24	2020/3/23	25	+
14	Nasopharyngeal swabs	2020/2/11	2020/2/10	26	+++
15	Nasopharyngeal swabs	2020/2/14	2020/1/8	30	±
16	Nasopharyngeal swabs	2020/2/9	2020/1/8	27	±
17	Nasopharyngeal swabs	2020/2/10	2020/1/27	26	+
18	Nasopharyngeal swabs	2020/1/28	2020/1/18	27	+
19	Nasopharyngeal swabs	2020/3/10	2020/3/7	26	+
20	Nasopharyngeal swabs	2020/2/14	2020/1/27	26	+
21	Nasopharyngeal swabs	2020/2/13	2020/2/10	27	+
22	Nasopharyngeal swabs	2020/1/31	2020/1/23	28	±
23	Nasopharyngeal swabs	2020/2/5	2020/1/28	27.5	+
24	Nasopharyngeal swabs	2020/2/23	2020/2/17	26	++
25	Nasopharyngeal swabs	2020/2/16	2020/2/10	25.5	+++
26	Nasopharyngeal swabs	2020/2/15	2020/2/5	29	+
27	Nasopharyngeal swabs	2020/1/25	2020/1/21	27	+

28	Nasopharyngeal swabs	2020/3/22	2020/3/17	25	++
29	Nasopharyngeal swabs	2020/9/27	2020/9/21	27	+
30	Nasopharyngeal swabs	2020/11/2	2020/10/29	27	±
31	Nasopharyngeal swabs	2020/9/27	2020/9/22	27	+
32	Nasopharyngeal swabs	2020/10/20	2020/10/16	25	+
33	Nasopharyngeal swabs	2020/11/2	2020/10/31	27	+
34	Nasopharyngeal swabs	2020/9/8	2020/9/2	25	++
35	Nasopharyngeal swabs	2020/8/16	2020/8/11	23	++
36	Nasopharyngeal swabs	2020/4/26	2020/4/19	25	++
37	Nasopharyngeal swabs	2020/4/23	2020/4/10	27	-
38	Nasopharyngeal swabs	2020/5/31	2020/5/24	27	±
39	Nasopharyngeal swabs	2020/5/6	2020/4/30	28	±
40	Nasopharyngeal swabs	2020/5/24	2020/5/20	27.5	+
41	Nasopharyngeal swabs	2020/5/13	2020/5/11	27	+
42	Nasopharyngeal swabs	2020/7/22	2020/7/17	25	++
43	Nasopharyngeal swabs	2020/7/9	2020/7/3	27	+

44	Nasopharyngeal swabs	2020/7/29	2020/7/23	27	±
45	Nasopharyngeal swabs	2020/6/28	2020/6/25	23	++
46	Nasopharyngeal swabs	2020/8/13	2020/8/7	23	+++
47	Nasopharyngeal swabs	2020/8/22	2020/8/19	25	+
48	Nasopharyngeal swabs	2020/6/27	2020/6/23	26	+++
49	Nasopharyngeal swabs	2020/6/1	2020/5/24	30	±
50	Nasopharyngeal swabs	2020/6/6	2020/6/2	27	±
51	Nasopharyngeal swabs	2020/8/17	2020/8/11	26	+
52	Nasopharyngeal swabs	2020/9/2	2020/8/26	22	+++
53	Nasopharyngeal swabs	2020/4/23	2020/4/18	22	++++
54	Nasopharyngeal swabs	2020/5/27	2020/5/22	26	+
55	Nasopharyngeal swabs	2020/5/14	2020/5/9	27	+
56	Nasopharyngeal swabs	2020/3/28	2020/3/23	26	+
57	Nasopharyngeal swabs	2020/4/12	2020/4/7	26	+
58	Nasopharyngeal swabs	2020/6/7	2020/6/1	27	+
59	Nasopharyngeal swabs	2020/5/31	2020/5/24	28	±

60	Nasopharyngeal swabs	2020/3/22	2020/3/19	24.5	+++
61	Nasopharyngeal swabs	2020/2/18	2020/2/13	27	+
62	Nasopharyngeal swabs	2020/4/17	2020/4/11	28	±
63	Nasopharyngeal swabs	2020/8/28	2020/8/23	25	+
64	Nasopharyngeal swabs	2020/8/21	2020/8/17	26	++
65	Nasopharyngeal swabs	2020/8/9	2020/8/3	25.5	++
66	Nasopharyngeal swabs	2020/9/27	2020/9/21	23	+++
67	Nasopharyngeal swabs	2020/9/29	2020/9/22	25	+
68	Nasopharyngeal swabs	2020/8/26	2020/8/21	26	+++
69	Nasopharyngeal swabs	2020/7/21	2020/7/14	30	+
70	Nasopharyngeal swabs	2020/7/28	2020/7/23	27.5	+
71	Nasopharyngeal swabs	2020/5/31	2020/5/25	25	+
72	Nasopharyngeal swabs	2020/5/22	2020/5/16	29	+
73	Nasopharyngeal swabs	2020/6/27	2020/6/19	27	+
74	Nasopharyngeal swabs	2020/5/14	2020/5/7	25	+
75	Nasopharyngeal swabs	2020/4/30	2020/4/26	27	+

76	Nasopharyngeal swabs	2020/7/28	2020/7/20	25	++
77	Nasopharyngeal swabs	2020/7/19	2020/7/14	26	+
78	Nasopharyngeal swabs	2020/8/30	2020/8/23	26	+
79	Nasopharyngeal swabs	2020/5/25	2020/5/19	22	++++
80	Nasopharyngeal swabs	2020/4/29	2020/4/22	24.5	++
81	Nasopharyngeal swabs	2020/6/21	2020/6/20	23	+++
82	Nasopharyngeal swabs	2020/6/12	2020/6/7	23	++
83	Nasopharyngeal swabs	2020/9/18	2020/9/13	22	++++
84	Nasopharyngeal swabs	2020/8/20	2020/8/17	26	+
85	Nasopharyngeal swabs	2020/8/1	2020/7/27	27	+
86	Nasopharyngeal swabs	2020/9/19	2020/9/11	26	+
87	Nasopharyngeal swabs	2020/7/22	2020/7/13	26	+
88	Nasopharyngeal swabs	2020/5/16	2020/5/11	28	-
89	Nasopharyngeal swabs	2020/5/5	2020/4/28	27.5	+
90	Nasopharyngeal swabs	2020/6/14	2020/6/12	27	+
91	Nasopharyngeal swabs	2020/7/27	2020/7/20	25	++

92	Nasopharyngeal swabs	2020/7/13	2020/7/8	25	++
93	Nasopharyngeal swabs	2020/6/28	2020/6/25	27	+
94	Nasopharyngeal swabs	2020/8/9	/	27	±
95	Nasopharyngeal swabs	2020/8/16	/	23	++
96	Nasopharyngeal swabs	2020/6/22	/	23	+++
97	Nasopharyngeal swabs	2020/9/10	/	25	++
98	Nasopharyngeal swabs	2020/10/19	/	27	+
99	Nasopharyngeal swabs	2020/10/29	/	22	++++
100	Nasopharyngeal swabs	2020/11/16	/	23	++++
何少欣	Nasopharyngeal swabs	2020/12/4	/	/	-
乔雪丽	Nasopharyngeal swabs	2020/12/4	/	/	-
黄炯桥	Nasopharyngeal swabs	2020/12/4	/	/	-
高帅鹏	Nasopharyngeal swabs	2020/12/4	/	/	-
孙小凤	Nasopharyngeal swabs	2020/12/4	/	/	-
方桂军	Nasopharyngeal swabs	2020/12/4	/	/	-
殷涛	Nasopharyngeal swabs	2020/12/4	/	/	-

覃健婷	Nasopharyngeal swabs	2020/12/4	/	/	-
廖带娣	Nasopharyngeal swabs	2020/12/4	/	/	-
汪光海	Nasopharyngeal swabs	2020/12/4	/	/	-
唐铭健	Nasopharyngeal swabs	2020/12/4	/	/	-
王婷婷	Nasopharyngeal swabs	2020/12/4	/	/	-
赵亚蒙	Nasopharyngeal swabs	2020/12/4	/	/	-
王丽兰	Nasopharyngeal swabs	2020/12/4	/	/	-
陈炜光	Nasopharyngeal swabs	2020/12/4	/	/	-
陈莉莉	Nasopharyngeal swabs	2020/12/4	/	/	-
黎耶玲	Nasopharyngeal swabs	2020/12/4	/	/	-
李虹倩	Nasopharyngeal swabs	2020/12/4	/	/	-
韦忠律	Nasopharyngeal swabs	2020/12/4	/	/	-
曹焯	Nasopharyngeal swabs	2020/12/4	/	/	-
农艳宁	Nasopharyngeal swabs	2020/12/4	/	/	-
龙洁云	Nasopharyngeal swabs	2020/12/4	/	/	-
杨路君	Nasopharyngeal swabs	2020/12/4	/	/	-

罗艳化	Nasopharyngeal swabs	2020/12/4	/	/	-
林艳	Nasopharyngeal swabs	2020/12/4	/	/	-
凌远佳	Nasopharyngeal swabs	2020/12/4	/	/	-
黎一锋	Nasopharyngeal swabs	2020/12/4	/	/	-
黄怡欣	Nasopharyngeal swabs	2020/12/4	/	/	-
黎燕	Nasopharyngeal swabs	2020/12/4	/	/	-
韦俏婷	Nasopharyngeal swabs	2020/12/4	/	/	-
黄婉	Nasopharyngeal swabs	2020/12/4	/	/	-
梁梅清	Nasopharyngeal swabs	2020/12/4	/	/	-
叶晓锋	Nasopharyngeal swabs	2020/12/4	/	/	-
伍锦池	Nasopharyngeal swabs	2020/12/4	/	/	-
林伟鸿	Nasopharyngeal swabs	2020/12/4	/	/	-
薛辉辉	Nasopharyngeal swabs	2020/12/4	/	/	-
李松宴	Nasopharyngeal swabs	2020/12/4	/	/	-
柯强强	Nasopharyngeal swabs	2020/12/4	/	/	-
杨宏委	Nasopharyngeal swabs	2020/12/4	/	/	-

曾剑锋	Nasopharyngeal swabs	2020/12/4	/	/	-
吴华堂	Nasopharyngeal swabs	2020/12/4	/	/	-
罗强	Nasopharyngeal swabs	2020/12/4	/	/	-
龙正章	Nasopharyngeal swabs	2020/12/4	/	/	-
罗龙俊	Nasopharyngeal swabs	2020/12/4	/	/	-
苏俊铨	Nasopharyngeal swabs	2020/12/4	/	/	-
张程	Nasopharyngeal swabs	2020/12/4	/	/	-
曾小明	Nasopharyngeal swabs	2020/12/4	/	/	-
王伟	Nasopharyngeal swabs	2020/12/4	/	/	-
张志勇	Nasopharyngeal swabs	2020/12/4	/	/	-
叶锦文	Nasopharyngeal swabs	2020/12/4	/	/	-
宋军	Nasopharyngeal swabs	2020/12/4	/	/	-
潘锦肖	Nasopharyngeal swabs	2020/12/4	/	/	-
张伟帆	Nasopharyngeal swabs	2020/12/4	/	/	-
张鹏飞	Nasopharyngeal swabs	2020/12/4	/	/	-
王小龙	Nasopharyngeal swabs	2020/12/4	/	/	-

叶德军	Nasopharyngeal swabs	2020/12/4	/	/	-
钟林威	Nasopharyngeal swabs	2020/12/4	/	/	-
陈炳奇	Nasopharyngeal swabs	2020/12/4	/	/	-
陈泽坤	Nasopharyngeal swabs	2020/12/4	/	/	-
周木炉	Nasopharyngeal swabs	2020/12/4	/	/	-
何嘉磊	Nasopharyngeal swabs	2020/12/4	/	/	-
蓝青扬	Nasopharyngeal swabs	2020/12/4	/	/	-
洪乐	Nasopharyngeal swabs	2020/12/4	/	/	-
唐正	Nasopharyngeal swabs	2020/12/4	/	/	-
蒙家铭	Nasopharyngeal swabs	2020/12/4	/	/	-
杨子航	Nasopharyngeal swabs	2020/12/4	/	/	-
陈建奇	Nasopharyngeal swabs	2020/12/4	/	/	-
温兆权	Nasopharyngeal swabs	2020/12/4	/	/	-
陆顺东	Nasopharyngeal swabs	2020/12/4	/	/	-
翁石林	Nasopharyngeal swabs	2020/12/4	/	/	-
刘勇	Nasopharyngeal swabs	2020/12/4	/	/	-

张凯林	Nasopharyngeal swabs	2020/12/4	/	/	-
张淡林	Nasopharyngeal swabs	2020/12/4	/	/	-
余志刚	Nasopharyngeal swabs	2020/12/4	/	/	-
徐泳烙	Nasopharyngeal swabs	2020/12/4	/	/	-
魏荟迺	Nasopharyngeal swabs	2020/12/4	/	/	-
潘光贤	Nasopharyngeal swabs	2020/12/4	/	/	-
鲁存	Nasopharyngeal swabs	2020/12/4	/	/	-
梁辉	Nasopharyngeal swabs	2020/12/4	/	/	-
吴焕森	Nasopharyngeal swabs	2020/12/4	/	/	-
罗轩	Nasopharyngeal swabs	2020/12/4	/	/	-
杨广浩	Nasopharyngeal swabs	2020/12/4	/	/	-
马国强	Nasopharyngeal swabs	2020/12/4	/	/	-
陈佳明	Nasopharyngeal swabs	2020/12/4	/	/	-
黄志聪	Nasopharyngeal swabs	2020/12/4	/	/	-
蒋鹏	Nasopharyngeal swabs	2020/12/4	/	/	-
裴俊良	Nasopharyngeal swabs	2020/12/4	/	/	-

莫海湘	Nasopharyngeal swabs	2020/12/4	/	/	-
廖嘉鑫	Nasopharyngeal swabs	2020/12/4	/	/	-
李 军	Nasopharyngeal swabs	2020/12/4	/	/	-
林观生	Nasopharyngeal swabs	2020/12/4	/	/	-
曹立华	Nasopharyngeal swabs	2020/12/4	/	/	-
王紫伦	Nasopharyngeal swabs	2020/12/4	/	/	-
林润喜	Nasopharyngeal swabs	2020/12/4	/	/	-
霍贝贝	Nasopharyngeal swabs	2020/12/4	/	/	-
李维煜	Nasopharyngeal swabs	2020/12/4	/	/	-
曾燕康	Nasopharyngeal swabs	2020/12/4	/	/	-
吴瀚选	Nasopharyngeal swabs	2020/12/4	/	/	-
廖慧威	Nasopharyngeal swabs	2020/12/4	/	/	-
罗培良	Nasopharyngeal swabs	2020/12/4	/	/	-
李堪贵	Nasopharyngeal swabs	2020/12/4	/	/	-
李昌明	Nasopharyngeal swabs	2020/12/4	/	/	-
邹进平	Nasopharyngeal swabs	2020/12/4	/	/	-

孙华檬	Nasopharyngeal swabs	2020/12/4	/	/	-
严立海	Nasopharyngeal swabs	2020/12/4	/	/	-
李于康	Nasopharyngeal swabs	2020/12/4	/	/	-
郑嘉莹	Nasopharyngeal swabs	2020/12/4	/	/	-
刘杰	Nasopharyngeal swabs	2020/12/4	/	/	-
张靖	Nasopharyngeal swabs	2020/12/4	/	/	-
黄荣发	Nasopharyngeal swabs	2020/12/4	/	/	-
许迪朗	Nasopharyngeal swabs	2020/12/4	/	/	-
叶伟军	Nasopharyngeal swabs	2020/12/4	/	/	-
曾启山	Nasopharyngeal swabs	2020/12/4	/	/	-
刘毅	Nasopharyngeal swabs	2020/12/4	/	/	-
温立基	Nasopharyngeal swabs	2020/12/4	/	/	-
宁磊	Nasopharyngeal swabs	2020/12/4	/	/	-
李江	Nasopharyngeal swabs	2020/12/4	/	/	-
阳基友	Nasopharyngeal swabs	2020/12/4	/	/	-
曾颖琪	Nasopharyngeal swabs	2020/12/4	/	/	-

杨永郁	Nasopharyngeal swabs	2020/12/4	/	/	-
韦炳锋	Nasopharyngeal swabs	2020/12/4	/	/	-
王凤琼	Nasopharyngeal swabs	2020/12/4	/	/	-
何玲东	Nasopharyngeal swabs	2020/12/4	/	/	-
欧阳植林	Nasopharyngeal swabs	2020/12/4	/	/	-
黄开金	Nasopharyngeal swabs	2020/12/4	/	/	-
伍菊青	Nasopharyngeal swabs	2020/12/4	/	/	-
肖兴连	Nasopharyngeal swabs	2020/12/4	/	/	-
黄辉文	Nasopharyngeal swabs	2020/12/4	/	/	-
陆玲	Nasopharyngeal swabs	2020/12/4	/	/	-
张海生	Nasopharyngeal swabs	2020/12/4	/	/	-
刘立芳	Nasopharyngeal swabs	2020/12/4	/	/	-
周爱群	Nasopharyngeal swabs	2020/12/4	/	/	-
张兰兰	Nasopharyngeal swabs	2020/12/4	/	/	-
王奕华	Nasopharyngeal swabs	2020/12/4	/	/	-
李雄艳	Nasopharyngeal swabs	2020/12/4	/	/	-

周佳惠	Nasopharyngeal swabs	2020/12/4	/	/	-
李羽菊	Nasopharyngeal swabs	2020/12/4	/	/	-
钟活其	Nasopharyngeal swabs	2020/12/4	/	/	-
余万河	Nasopharyngeal swabs	2020/12/4	/	/	-
吴丽敏	Nasopharyngeal swabs	2020/12/4	/	/	-
戴雪梅	Nasopharyngeal swabs	2020/12/4	/	/	-
杜姣姣	Nasopharyngeal swabs	2020/12/4	/	/	-
唐翠华	Nasopharyngeal swabs	2020/12/4	/	/	-
陈苇	Nasopharyngeal swabs	2020/12/4	/	/	-
黄陈潇	Nasopharyngeal swabs	2020/12/4	/	/	-
王明富	Nasopharyngeal swabs	2020/12/4	/	/	-
缪斯琼	Nasopharyngeal swabs	2020/12/4	/	/	-
卢玲玲	Nasopharyngeal swabs	2020/12/4	/	/	-
杨聪	Nasopharyngeal swabs	2020/12/4	/	/	-
韩帅杰	Nasopharyngeal swabs	2020/12/4	/	/	-
杨婉玲	Nasopharyngeal swabs	2020/12/4	/	/	-

李满红	Nasopharyngeal swabs	2020/12/4	/	/	-
关文琼	Nasopharyngeal swabs	2020/12/4	/	/	-
罗娅	Nasopharyngeal swabs	2020/12/4	/	/	-
李嘉灵	Nasopharyngeal swabs	2020/12/4	/	/	-
杨嘉满	Nasopharyngeal swabs	2020/12/4	/	/	-
陈润烈	Nasopharyngeal swabs	2020/12/4	/	/	-
李阳	Nasopharyngeal swabs	2020/12/4	/	/	-
罗晓凤	Nasopharyngeal swabs	2020/12/4	/	/	-
郭景辉	Nasopharyngeal swabs	2020/12/4	/	/	-
莫嘉弘	Nasopharyngeal swabs	2020/12/4	/	/	-
林静婷	Nasopharyngeal swabs	2020/12/4	/	/	-
林金帆	Nasopharyngeal swabs	2020/12/4	/	/	-
黄春智	Nasopharyngeal swabs	2020/12/4	/	/	-
蒋平	Nasopharyngeal swabs	2020/12/4	/	/	-
李娜	Nasopharyngeal swabs	2020/12/4	/	/	-
冯彭玲	Nasopharyngeal swabs	2020/12/4	/	/	-

倪泽慧	Nasopharyngeal swabs	2020/12/4	/	/	-
周玲	Nasopharyngeal swabs	2020/12/4	/	/	-
袁国飞	Nasopharyngeal swabs	2020/12/4	/	/	-
袁永佳	Nasopharyngeal swabs	2020/12/4	/	/	-
赵燕	Nasopharyngeal swabs	2020/12/4	/	/	-
贺奕恩	Nasopharyngeal swabs	2020/12/4	/	/	-
陆采红	Nasopharyngeal swabs	2020/12/4	/	/	-
严思萍	Nasopharyngeal swabs	2020/12/4	/	/	-
白思雨	Nasopharyngeal swabs	2020/12/4	/	/	-
余子灏	Nasopharyngeal swabs	2020/12/4	/	/	-
罗晓棠	Nasopharyngeal swabs	2020/12/4	/	/	-
赖文升	Nasopharyngeal swabs	2020/12/4	/	/	-
盘献芬	Nasopharyngeal swabs	2020/12/4	/	/	-
张杰	Nasopharyngeal swabs	2020/12/4	/	/	-
张俊	Nasopharyngeal swabs	2020/12/4	/	/	-
庞惠洁	Nasopharyngeal swabs	2020/12/4	/	/	-

韦凤霞	Nasopharyngeal swabs	2020/12/4	/	/	-
罗晓娇	Nasopharyngeal swabs	2020/12/4	/	/	-
林源	Nasopharyngeal swabs	2020/12/4	/	/	-
尹翔	Nasopharyngeal swabs	2020/12/4	/	/	-
唐西江	Nasopharyngeal swabs	2020/12/4	/	/	-
邱蕊萍	Nasopharyngeal swabs	2020/12/4	/	/	-
吕小玲	Nasopharyngeal swabs	2020/12/4	/	/	-
王璇	Nasopharyngeal swabs	2020/12/4	/	/	-
赖浩鑫	Nasopharyngeal swabs	2020/12/4	/	/	-
詹智豪	Nasopharyngeal swabs	2020/12/4	/	/	-
冉敏	Nasopharyngeal swabs	2020/12/4	/	/	-
郭嘉明	Nasopharyngeal swabs	2020/12/4	/	/	-
李广	Nasopharyngeal swabs	2020/12/4	/	/	-
严伟华	Nasopharyngeal swabs	2020/12/4	/	/	-
仲浩	Nasopharyngeal swabs	2020/12/4	/	/	-
林静怡	Nasopharyngeal swabs	2020/12/4	/	/	-



陈佩宜	Nasopharyngeal swabs	2020/12/4	/	/	-
刘韩池	Nasopharyngeal swabs	2020/12/4	/	/	-
唐莎	Nasopharyngeal swabs	2020/12/4	/	/	-
李文杰	Nasopharyngeal swabs	2020/12/4	/	/	-
史兴远	Nasopharyngeal swabs	2020/12/4	/	/	-
唐雅娴	Nasopharyngeal swabs	2020/12/4	/	/	-
黄宇豪	Nasopharyngeal swabs	2020/12/4	/	/	-
伍仪	Nasopharyngeal swabs	2020/12/4	/	/	-
曹宇	Nasopharyngeal swabs	2020/12/4	/	/	-
田宝磊	Nasopharyngeal swabs	2020/12/4	/	/	-
叶主衡	Nasopharyngeal swabs	2020/12/4	/	/	-

样品编号	Ct 值	名称	交叉反应
J01	27	冠状病毒 HKU1	-
J02	24	冠状病毒 OC43	-
J03	28	冠状病毒 NL63	-
J04	31	冠状病毒 229E	-
J05	27	HPIV-1	-
J06	30	HPIV-2	-
J07	34	HPIV-3	-
J08	24	HPIV-4	-
J09	31	HRV	-
J10	25	HAV	-
J11	29	HRSV	-
J12	20	肺炎支原体	-

